



**PREVENTION** is better than cure, what if you could avoid ACD\* in the first place?

\*Allergic Contact Dermatitis (ACD) Type IV: Delayed Contact Dermatitis



Itchiness



Redness



Peeling Skin

## SafeTouch Advanced Pure

Accelerator Free Medical Examination Gloves

The **Medicom SafeTouch Advanced Pure** is the ideal glove solution for sensitive hands. These gloves are manufactured without the use of the chemical accelerators resulting in a cleaner, more skin-friendly product. It provides the best possible solution when you need the barrier protection of a glove combined with healthy skin protection for sensitive hands.



## SafeTouch Advanced Pure | Key Features

 <b>Allergy Care</b> ASTM D7558-09 ISO 10993-10	 <b>Medical Grade</b> ASTM D6319	 <b>Protection</b> EN 420 EN ISO 374
---	---	--

 Aged Care	 Beauty
 Dental	 Food & Agriculture
 Industrial	 Life Science
 Medical	

- Accelerator Free 
- Medical Grade 
- Chemo Drugs Tested 
- Food Grade 
- Fingertip Textured 

- DermaSense Technology  
Skin Friend product 
- Ambidextrous 
- Beaded Cuff 
- Powder Free 
- Latex Free 

### Product Code

<b>UN131275A</b>	XS
<b>UN131275B</b>	S
<b>UN131275C</b>	M
<b>UN131275D</b>	L
<b>UN131275E</b>	XL



### Packaging

10 boxes x 200 pcs

### Length

Min. 240mm

### Colour



**Proud Leader**  
in Infection Control  
Since 1988

### Contact Us:

Tel: 1300 818 168  
Email: [cs@medicom-australia.com](mailto:cs@medicom-australia.com)  
Visit: [au.medicom.com](http://au.medicom.com)